Lava™

Precision Solutions
精密修复的最佳系统

www.dentalab.com

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3M

Lava简介

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About Lava

Lava is the brand name of 3M digital dentistry. It provides dentists a total precision solution by a digital workflow. Lava includes Lava C.O.S (Chairside Oral Scanner), Lava CAD/CAM system, Lava Connect network and Lava machinable materials. Currently, Lava CAD/CAM system, Lava Connect network and Lava Zirconia have been launched in China market. In near future, a series of machinable materials will enter China, i.e. ultra translucent zirconia, glass ceramic and nano resin composite.

Since its introduction in 2001, Lava Zirconia has become a huge success story. Millions of restorations have been produced and Lava System are running in 40 countries around the globe.

Lava Zirconia is one of the most thoroughly tested CAD/CAM materials on the market with over 100 studies published by researchers world wide. Clinical excellence is proven by more than 10 clinical studies with more than 1,500 restorations placed and followed up to eight years. These studies validate why Lava is known for high-strength with outstanding marginal fit and excellent aesthetics.

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Are all zirconias the same?

Ceramic restorations in the posterior region were once limited to single units. CAD/CAM technology made it possible to prepare restorations out of high-strength ceramics. Now, with the introduction of zirconia as a dental material, clinicians can place multi-unit restorations in both the anterior and posterior regions. This is due in part to the high flexural strength and fracture toughness of zirconia. Zirconia materials from different manufacturers, however, may be processed differently and have varying levels of stability. Not all manufacturers, for example, have completed adequate in vitro and in vivo clinical studies. Final restoration quality is directly dependent on careful and accurate control of the manufacturing process and thorough testing to substantiate material reliability.

What can be different?


1. Processing parameters for pre-sintered zirconia affect performance attributes.
2. Differences in the zirconia powder affect the strength/long-term stability and translucency of the restoration.
3. The pressing condition and pressing method affect the marginal fit, strength and translucency of the restoration.
4. Pre-sintering conditions affect the strength of the pre-sintered material and its millability.
5. Coloring of the zirconia can affect the marginal fit, strength and translucency of the material.

Under multiple options for you, dentists should ask your partner lab for the brand, source, license; in vitro and in vivo research data and warranty.

Although the zirconia ceramic is chemically similar, once processed, it can exhibit different mechanical, optical characteristics and longevity.

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
Advantages of Zirconia and indication

For many years now, we have only had two basic options when it came to restorations--

- 1) Something strong but with esthetic limitations-- All metal or Porcelain fused to metal
- 2) Or something life-like and esthetic but very brittle

With the introduction of zirconia, we now are approaching "the best of both worlds" – strength and beauty. In the environment of oral cavity with complex stress and temperature alternation, zirconia has longer longevity as compared to other all-ceramics.

Therefore, zirconia has wider indication from single to multiunits and custom-made abutment.

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Lava Zirconia

As mentioned above, zirconia is not alike. Lava zirconia is manufactured from strict raw materials purification to final computerized sintering through a highly stringent control in process. Lava offers dentists restorations with incomparable high quality as followings,

1. High strength and long-term stability under both flexural loading and aging condition.


2. Perfect marginal fit

Every Lava Zirconia blank has a unique shrinkage factor coded into the blank to ensure precision process consistency.

3. Excellent translucency and patented shading system guarantee esthetics.

4. Lava™ Authentication Program.

5. The study with the largest case number and longest restoration time was published by THE DENTAL ADVISOR™. Lava was assessed to "perform exceptionally well".

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Is zirconia safe?

There are no known allergic reactions of medical device made of zirconia. Furthermore, zirconia reveals absolutely no potential for cytotoxicity.

Some patients may worry about the radioactivity of zirconia restorations. When the raw mineral materials are processed, the strict procedure will be implemented to purify them to remove off all impurities. All the zirconia for medical usage need fulfill the requirements of ISO standard, as well as other all-ceramics. It has been proved purified zirconia powder is absoluty safe. We can learnt that how important the quality, source and process of manufacturing are for your product selection.

Claims which maintain that zirconia is metal, are false. Such misinterpretations are probably because that people sometimes is confused by the terms zirconium or zircon and zirconia.

In conclusion, zirconia is safe to be used as dental ceramics. Lava Zirconia Frame is 'Made in Germany'. To be precise in the 3M ESPE facility in Seefeld, Bavaria. Here the material were optimized over all steps to achieve a perfect outcome, including safety.

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

Chipping of zirconia supported restoration

Recent studies demonstrated that there was no statistical difference of chipping between zirconia and metal supported restorations. Chipping in zirconia happened often internal structure of veneerings, but the interface. The major reasons includes :

- (1) improper coping design, thicker porcelain;
- (2) Inappropriate fusion control.
- (3) For dentists, no or incorrect further treatment after grinding will lead to chipping.

3 step polishing or glaze is mandatory.

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Is grinding detrimental to crown?

As mentioned above, adjustments resulted in surface flaws which can serve as crack initiation sites. If possible - Don't do it!

If you have to do adjustments:

- Take a very fine diamond bur (red ring or finer)
- Use copious amounts of water
- Polish, Polish, Polish or glazed

If you have to do large corrections on a regular basis

- Check your impressing, bite taking and cementation techniques
- Discuss with your lab
- Do a try-in and corrections BEFORE the final glaze firing

Once coping or framework is sintered and if any grinding is necessary, use a new fine diamond (e.g. Brasseler #8881) under copious amounts of water just for very limited area. Please contact your lab partner, if you have to grind zirconia coping and consider the possibility of fracture.

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Will all zirconia cause over abrasion of antagonist?

For the patients with little inter-occlusal space or bruxism, dentists can choose all zirconia restorations.

Summary of the actual literature:

- To minimize antagonist wear, the surface of the restorative material has to be fine-grained and smooth.
- If monolithic Zirconia is used, the surface must always be polished or glazed when the antagonist is natural tooth.
- Less enamel wear should be expected on polished surfaces due to the lower friction coefficient.
- Surface roughness of the antagonistic ceramic influenced the wear of human enamel

Recommendation of polisher brands: EVE, DFS, Komet

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Cementation of zirconia

Lava strength allows conventional cementation for crowns and bridges

You can make single posterior crowns out of feldspathic glass. However, these materials need to be adhesively bonded to the tooth to stabilize the restoration. Adhesive cementation can be challenging e.g. in cases with sub gingival margins, where excess removal is cumbersome.

RelyX Unicem contains phosphoric acid functional groups which bond to zirconia surface. It is highly recommended.

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Due to its excellent mechanical and optical characteristics, 3M™ ESPE™ Lava™ Zirconia frameworks can be used for a wide range of indications.

图1: 单冠
Fig. 1: Single Crowns

图2: 联冠
Fig. 2: Splinted Crowns

图3: 3单位桥
Fig. 3: 3-unit Bridges


图4: 4单位桥
Fig. 4: 4-unit Bridges


图5: 5单位桥
Fig. 5: 5-unit Bridges


图6: 6单位桥
Fig. 6: 6-unit Bridges

图7: 有曲度的和长跨度的桥(最长达48mm)
Fig. 7: Curved and Long-span Bridges

图8: 悬臂桥(磨牙症患者除外)
Fig. 8: Cantilever Bridges
(excluded for patients with bruxism)

图9: 嵌体桥和高嵌体桥(磨牙症患者除外)
Fig. 9: Inlay Bridges and Onlay Bridges
(excluded for patients with bruxism)

图10: 前牙粘结性桥(磨牙症患者除外)
Fig. 10: Anterior Adhesive Bridges
(excluded for patients with bruxism)

图11: 套筒冠
Fig. 11: Primary Crowns/Telescopes

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
Ideal Preparation: Shoulder or Chamfer?

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
A reduction of the tooth structure based on the dimensions indicated below is sufficient. We recommend a preparation matrix of the initial clinical situation in order to check the progress of the tooth preparation. Ideally, the preparation includes a circumferential shoulder or chamfer with a horizontal angle of at least 5°. The vertical preparation angle should be at least 4°. The inside angle of the shoulder preparation must be given a rounded contour. All occlusal and incisal edges should also be rounded. The marginal edge of the preparation needs to be continuous and clearly visible. A bevel should be avoided.

Ideal

前牙建议备牙量
Recommended preparation for anterior teeth.

后牙建议备牙量
Recommended preparation for posterior teeth.

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

沟槽式制备：边缘不能被清晰地检出。
Gutter Preparation: Margin cannot be detected clearly.

90° 肩台：边缘不能被清晰地检出。
90° Shoulder: Margin cannot be detected clearly.


一定要避免倒凹。
Undercuts must be avoided.

平行壁：原则上，平行壁制备是可行的。然而，在这种情况下粘结间隙是切割不出来的。这可能会在很大程度上影响就位。
Parallel wall preparations result in a cement gap that cannot be milled. This may significantly affect the fit. A 5° taper should be used.

一定要避免尖锐的切缘-殆面边缘。
弧形半径应>0.4mm。
Sharp incisal-occlusal edges must be avoided. The rounding radius should be > 0.4 mm.

由于就位道的限制，不能制作有倾斜牙的桥。
Bridges with inclined teeth cannot be achieved due to the restricted path of insertion.